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USER'S MANUAL

READ AND UNDERSTAND THIS MANUAL
BEFORE OPERATION. SAVE THIS MANUAL.

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SHIMADZU CORPORATION
ANALYTICAL INSTRUMENTS DIVISION
KYOTO, JAPAN

Copyright © Shimadzu Corporation 1992. All rights are reserved, including those to reproduce this publication or parts thereof in any form without notice and does not represent a commitment on the part of the vendor.

Any errors or omissions which may have occurred in this publication despite the utmost care taken in its production will be corrected as soon as possible, but not necessarily immediately upon detection.

Shimadzu Corporation will continue to supply original Shimadzu maintenance parts associated with a given product for a period up to 10 years from the cessation of production of that product. Please be informed that after this period of time, Shimadzu Corporation cannot guarantee supply of such original maintenance parts. However, following discontinuation of a product, the period of availability of maintenance parts which have been produced on a sub-contract basis is up to the discretion of the concerned sub-contracting company.

Precautions for Safe Operation

This drain tube kit is an accessories of the Shimadzu high performance liquid chromatograph series LC-10A.

To operate the kit safely, strictly observe the following precautions.

1. Use the unit only for the purpose mentioned above.
2. Follow the procedures described in the user's manual.
3. Observe all warnings and cautions.
4. Do not disassemble or modify the unit except that the procedure is described in the user's manual.
5. Repair and Maintenance of the Unit.

Normal maintenance of this unit can be performed without removing the cover. Do not remove the cover at the normal maintenance.

Contact our sales office or agent if repair requires to remove the cover of the main body.

6. Installation of this unit may be dangerous. Only qualified Shimadzu service personnel may install this unit; the user is prohibited from attempting installation. The pink-color pages and the shaded parts are intended solely for the use of Shimadzu service personnel.
7. This system has the mark  at the portion where special caution is required.

When working at the portion indicated with this mark, read the instruction manual and exercise full caution.

MEANING OF CAUTION SIGNS

CAUTION	indicates a potentially hazardous situation which, if not avoided, may result in minor injury, or may result in machine damage.
----------------	---

Precaution on installation site and handling of the unit

High performance liquid chromatography, which often requires large volumes of organic solvent, requires proper care in handling and appropriate installation site.

Be sure to pay due attention to the precautions noted throughout manual, in addition to those below.

1. Ventilation

Solvents used in high performance liquid chromatography are often inflammable and toxic, therefore the room in which they are used must be well ventilated.

2. Fire

Never use an open flame in the vicinity of a high performance liquid chromatograph. Do not install other equipment that may generate sparks in the same room.

Install a fire extinguisher at a conspicuous location nearby.

3. Protective Goggles

Wear protective goggles when using solvents.

4. Other Facilities Required

A sink should be located as close as possible to the unit so that toxic solvents splashed in the eyes or on the skin may be immediately rinsed off.

Precautions on Static Electricity

Liquid chromatography using flammable organic solvents as mobile phase requires proper care against fire, explosion, etc. Particularly, among various possible accidents, those caused by static electricity are difficult to anticipate, and tend to occur only with unexpected conditions which often make countermeasures insufficient.

At a site where preparative liquid chromatography is practiced, a large amount of flammable substances may be used. Therefore, once an accident happens, it could lead to tremendous damage. The mechanism of accident caused by static electrical discharge and preventive measures are described below. Take due care in safety measures in handling of equipment.

1. Mechanism of Static Electrical Discharge Accident (Example)

Accidents caused by static electricity take place through the following processes.

Occurrence of Static Electricity

When liquid is fed at high speed through a small-diameter tube like the pipe of a liquid chromatograph, static electrical charge occurs by friction between solid and liquid as shown in Fig.12.1.


Fig.12.1 Occurrence of Static Electricity by Friction between Solid and Liquid

Charging and storage of static electricity

Energy release by discharge

Ignition of combustible substances

When the charged liquid is collected in an insulated vessel, the static charge accumulates gradually, and the voltage can easily reach a few kilovolts.

If some other conductive object is brought near the vessel, electricity is discharged at a certain distance from the vessel releasing heat energy.

If flammable gas of sufficient concentration exists nearby, ignition is caused by this energy.


Fig.12.2 Conditions for Accidents

2. Preventive Measures against Accidents

The principal preventive measure is the prevention of "charging and storage of static electricity" among those items shown in "Occurrence mechanism of static electrical discharge accidents." The preventive measures are shown below. It is recommended to exercise two or more measures simultaneously.

※Particularly when a large quantity of flammable solvent is held in a large vessel, be sure to observe the preventive measures 1, 2, and 3.

Preventive measure 1.

Use metallic (conductive) waste liquid vessel which is well grounded. This releases the charge of the waste liquid and vessel to ground.

The following items are available.

- (1) Grounding wire with clip P/N 228-21353-91
- (2) Metallic 18 liter can P/N 038-00044
- (3) Metallic 4 liter can P/N 038-00043-01

※Be sure to ground the vessel properly. Disconnecting of grounding wire or poor grounding defeat the purpose of using a metallic vessel.

※There are some metallic cans which have no conductivity due to an oxidized coating or lacquer on their surface. Be sure to confirm the grounding of vessels by a tester before application.

※When a liquid with almost no conductivity (of 10^{-10} S/m or less) is discharged into the vessel, it is necessary to mix it with another liquid with some conductivity. (The other liquid can be placed in the vessel in advance.)

Preventive measure 2.

Minimize the clearance of both inlet and outlet of vessel to prevent flame from entering the vessel.

(1) Cap with three holes for 18 liter and 4 liter cans (P/N 228-21354-91) is available.


Fig.12.3 Anti-Static Electricity Measures for Vessel

Preventive measure 3.

Do not approach the vessel with charged objects including the human body.

Charging prevention measures for human body.

- a) Prevention of charging of shoes and clothes
- b) Grounding of human body
- c) Make working floor conductive

Suitable products to be used for those measures a), b), and c) are available on the market.

※When persons who use no charge prevention measures approach dangerous sections, they have to be grounded beforehand. (For example, they should contact grounded metal by hand.)

Preventive measure 4.

Use pipes with inner diameter of 2mm or more for waste liquid line for large flow rates.

※Inclusion of bubbles in the tube may increase the amount of charging by ten times.

Check that there is no inclusion of air via tube joints.

Preventive measure 5.

When it is impossible to use a conductive vessel, use caution in the following points.

- a) Set the vessel so that the pipe outlet will be placed below the liquid level in the vessel. Or, dip a grounded metal (ex. pipe connected to the main body of device) in the liquid.

※This method is not effective for liquid with small conductivity (10^{-8} s/m or less).

- b) Use a vessel of the smallest possible capacity to minimize the damage by fire if it should occur.
- c) Prevent the room from being dry. Humidity of 65% or more has charge prevention effects.

1. Outline

This drain tube kit is connected to each component of the Shimadzu high performance liquid chromatograph series LC-10A. If the liquid leaked inside the system, this kit introduces the leaked liquid to the drain bottle so that the liquid may not splash onto the other device or may not collect on the desk.


Caution:

Each component of the LC-10A series has structure to drain the liquid from the connecting port for the leaked liquid exhaust tube at the lower portion of the right side of the system.
When installing the system, arrange plumbing for this drain tube kit.

2. Parts list

The accessories consist of the following parts. Check that all parts shown below are packaged.

No.	Part name	Part No.	Q'ty
1	Drain tube	228-25162-03	1
2	L-type drain tube	228-28094	1
3	T-type joint	228-28162	1
4	Straight type joint	228-28163	1
5	Pipe clip	037-60177-05	1


Caution:

When plumbing and carrying out maintenance, be sure to use the parts described on this page. Use of other parts may not assure normal function of the system.

3. Connection method


- (1)Connect the L-type drain tube to the connecting port for leaked liquid exhaust tube of each component.
- (2)Connect the device at the botton to the drain tube via the straight type joint, then discharge the liquid to the drain bottle.
- (3)Connect the components located at the second step from the bottom to the drain tube coming from the component stacked above it via the T-type joint and introduce the drain tube to the drain bottle. Connect the components above the third step

from the bottom to the upper component and the lower component via the T-type joint. When the drain tube is too long for connection, cut it at a proper length. Additionally, attach the L-type drain tube downward from the horizontal level and fix it to the front of the system with the **pipe clip** (accessory) if necessary. Place the drain bottle at the position lower than the device at the bottom.


Caution:

In order to prevent backflow of the liquid, fix the tube at the position lower than the drain port of the devices so that the liquid may flow smoothly.

Caution:

If the flow line of leaked liquid exhaust tube is placed above the drain port, the liquid cannot be drained. Be sure to place the waste water bottle at the position lower than the system (like on the floor) and make plumbing of the drain tube downward. After plumbing, flow the water to the receiving pan to check that the waste liquid does not flow back and is introduced to the drain bottle.

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